

# EntericBio<sup>®</sup>

## Revolutionising Molecular Enteric Testing



- Test Direct From Sample, No DNA/RNA Extraction Or Purification
- Freeze Dried, Room Temp Stable & Ready To Use Assays
- Automated Setup With Pierceable Cap Technology
- Easy To Use
- Minimal Hands On Time
- High Throughput

The system can process 1-48 samples per run with up to 276 patient samples in under 8 hours with no DNA/RNA extraction, automated setup and lyophilised ready to use assays.

fast, reliable & easy to use

[entericbio.com](http://entericbio.com)



Add Samples



Heat Treat



Automated Setup



Amplify & Detect



EntericBio®

	Combined		Bacteria				Parasite	Viral			CPE
	Dx	Gastro Panel 2	Gastro Panel 1	Bacterial Panel 3	Bacterial Panel 2	C.difficile	Parasite Panel 1	Viral Panel 3	Viral Panel 1	Viral Panel 2	CPE Screen
Salmonella	✓	✓	✓	✓							
Shigella spp. / EIEC*	✓	✓	✓	✓							
Campylobacter	✓	✓	✓	✓							
STEC* (stx 1/2)	✓	✓	✓	✓	✓						
Yersinia	✓			✓							
Vibrio	✓			✓							
Cryptosporidium	✓	✓					✓				
Giardia	✓	✓					✓				
Entamoeba histolytica	✓						✓				
C.difficile						✓					
Adenovirus								✓	✓		
Rotavirus								✓	✓		
Sapovirus								✓	✓		
Astrovirus								✓		✓	
Norovirus GI								✓		✓	
Norovirus GII								✓		✓	
OXA 48											✓
KPC											✓
NDM											✓
VIM											✓
IMP											✓
GES											✓

\*EIEC: Enteroinvasive E. Coli

\*STEC: Shiga toxin-producing E. Coli

EntericBio.com

Serosep Limited

Annacotty Business Park  
Annacotty, Limerick,  
Ireland

Telephone +353 (0)61 358 190  
Email info@serosep.com

Serosep UK Limited

14 The Brunel Centre  
Newton Road  
Crawley, West Sussex  
RH10 9TU, UK

Telephone +44 (0) 129 388 9215

No DNA / RNA Extraction • No Manual Pipetting • Rapid Results within 3 hrs • High throughput

